



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK
Tel: +44 (0)1792 686070
Email: info@enzymeresearch.co.uk
Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	BFIX 1005
Lot:	SAMPLE
Description:	Bovine Factor IX
Concentration:	1.44 mg/mL
Aliquots:	1 x 100 units
Volume/aliquot:	0.347 mL
Activity:	100.00 units
Total Protein:	0.50 mg (1 x 0.5 mg)
Buffer:	20 mM Tris-HCl / 0.1 M NaCl / 1mM Benzamidine / pH 7.4
Storage:	<-60 °C
Expiration Date:	** / ****

The Bovine Factor IXa beta was prepared from homogeneous Bovine Factor IX by activation with Bovine Factor XIa. This Bovine Factor XIa was removed after activation. Complete activation was observed on 10 % SDS-PAGE gels.

All bovine material is from U.S. cows <= 30 months of age.

For Research Use Only.