



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK
Tel: +44 (0)1792 686070
Email: info@enzymeresearch.co.uk
Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	BFX 1006
Lot:	SAMPLE
Description:	Bovine Factor X
Concentration:	1.77 mg/mL (before lyophilization)
Aliquots:	1 x 1.00 mg
Volume/aliquot:	0.565 mL
Activity:	N/A
Total Protein:	1.00 mg (1 x 1.0 mg)
Buffer:	20 mM Tris-HCl / 0.1 M NaCl / 1 mM Benzamidine / pH 7.4
Storage:	2-8 °C
Expiration Date:	**/****

The Bovine Factor X is >95% pure on SDS-PAGE and shows total reduction upon incubation with 2-mercaptoethanol.

All bovine material is from U.S. cows <= 30 months of age.

We suggest hydrating the protein with a good quality water to the original volume. After hydration aliquot into a useful (one time use) size and freeze at <-60°C.

For Research Use Only.