

Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK

Tel: +44 (0)1792 686070

Email: info@enzymeresearch.co.uk

Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	BPC 1003
Lot:	SAMPLE
Description:	Bovine Protein C
Concentration:	2.34 mg/mL (before lyophilisation)
Aliquots:	1 x 0.10 mg
Volume/aliquot:	0.043 mL
Activity:	N/A
Total Protein:	0.10 mg (1 x 0.1 mg)
Buffer:	20 mM Tris-HCl / 0.1 M NaCl / 1mM Benzamidine / pH 7.4
Storage:	2-8 °C
Expiration Date:	**/****

The Bovine Protein C is a >95% on SDS-PAGE gels and shows total reduction upon incubation with 2-mercaptoethanol. All bovine material is from U.S. cows <= 30 months of age.

We suggest hydrating the protein to the original volume with water. After hydration aliquot into a useful (one time use) size and freeze at <-60 °C

For Research Use Only.