



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK
Tel: +44 (0)1792 686070
Email: info@enzymeresearch.co.uk
Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Catalogue No:	Hir
Lot:	SAMPLE
Description:	r-Hirudin (HV-1)
Concentration:	n/a
Aliquots:	1 x 0.1 mg
Volume/aliquot:	n/a
Activity:	11,700 ATU/mg
Total Protein:	0.10 mg
Buffer:	n/a
Storage:	4 °C
Expiration Date:	**/****

The recombinant hirudin is >99% pure as judged by reverse phase hplc. No additives are included (buffer, salts etc) in the product.

For Research Use Only.