

Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK

Tel: +44 (0)1792 686070

Email: info@enzymeresearch.co.uk

Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	HFVIIa
Lot:	SAMPLE
Description:	Human Factor VIIa
Concentration:	1.04 mg/mL (before lyophilization)
Aliquots:	1 x 0.10 mg
Volume/aliquot:	0.096 mL
Activity:	53279.00 PEU/mg
Total Protein:	0.10 mg (1 x 0.1 mg)
Buffer:	0.03 M Tris-HCl, 0.15 M NaCl, 0.1% PEG 8000, 4% Ficoll Type 400, pH 7.4
Storage:	2-8°C
Expiration Date:	**/****

The Human Factor VIIa was prepared from homogeneous Human Factor VII using Human Factor XIIa. The Factor XIIa was removed using a Corn Trypsin Inhibitor column. The FVIIa is >95% pure by SDS-PAGE. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA and HBV DNA and Syphilis by FDA approved methods. Donors are screened for CJD.

We suggest hydrating the protein with a good quality water to the original volume. After hydration aliquot into a useful (one time use) size and freeze at <-60C.

For Research Use Only.