



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK
Tel: +44 (0)1792 686070
Email: info@enzymeresearch.co.uk
Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	HFXIIab
Lot:	SAMPLE
Description:	Human Factor XIIa beta
Concentration:	1.16 mg/mL (before lyophilization)
Aliquots:	15 x 0.10 mg
Volume/aliquot:	0.086 mL
Activity:	N/A
Total Protein:	1.50 mg (15 x 0.10 mg)
Buffer:	4 mM Sodium Acetate-HCl, 0.15 M NaCl, pH 5.3
Storage:	2-8 °C
Expiration Date:	**/****

The Human Factor XIIab was activated from homogeneous Human Factor XII using Dextran Sulfate. This Dextran Sulfate was removed after activation. Complete activation was observed on SDS-PAGE. The above protein was purified from Human plasma and found negative by FDA accepted methods for Anti-HIV 1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease). We suggest hydrating the protein to the original volume with water.

For Research Use Only.