



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK

Tel: +44 (0)1792 686070

Email: info@enzymeresearch.co.uk

Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	HPlasmin
Lot:	SAMPLE
Description:	Human Plasmin
Concentration:	1.32 mg/mL
Aliquots:	1 x 1.00 mg
Volume/aliquot:	0.758 mL
Activity:	7.40 U/mg (123.20 nkat/mg)
Total Protein:	1.00 mg
Buffer:	50 HEPES / 50mM Sodium Acetate / pH 8.5 / 50% Glycerol
Storage:	-20°C
Expiration Date:	**/****

The Human Plasmin was activated from homogeneous glu-plasminogen using urokinase. Complete activation was observed using the S-2251 chromogenic substrate. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA and HBV DNA and Syphilis by FDA approved methods. Donors are screened for CJD.