

Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK

Tel: +44 (0)1792 686070

Email: info@enzymeresearch.co.uk

Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Catalogue No:	2HK 1301
Lot:	SAMPLE
Description:	Two Chain HMW Kininogen
Concentration:	1.13 mg/mL (before lyophilization)
Aliquots:	4 x 1.00 mg
Volume/aliquot:	0.885 mL
Activity:	n/a
Total Protein:	4.00 mg (4 x 1.0 mg)
Buffer:	4 mM Sodium Acetate-HCl, 0.15 M NaCl, pH 5.3
Storage:	2-8 °C
Expiration Date:	**/****

The Two Chain High MW Kininogen is greater than 95% pure as judged by SDS-PAGE and shows complete reduction upon incubation with 2-mercaptoethanol. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

We suggest hydrating the protein to the original volume with water. Once hydrated, transfer the protein to a plastic vial. Unused portion may be aliquoted and frozen at <-60°C.

For Research Use Only.