



Enzyme
Research
Laboratories

Global Coag

Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK

Tel: +44 (0)1792 686070

Email: info@enzymeresearch.co.uk

Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	Mab HFXII
Lot:	SAMPLE
Description:	Monoclonal antibody to Human Factor XII
Concentration:	1 mg/mL (before lyophilization)
Aliquots:	1 x 0.50 mg
Volume/aliquot:	0.500 ml
Activity:	N/A
Total Protein:	0.50 mg
Buffer:	10 mM Sodium Phosphate / 0.15 M NaCl / 0.01% Sodium Azide / pH 7.4
Storage:	-20°C
Expiration Date:	**/****
Immunogen:	Human Factor XII
Clone ID:	GMA-140

Antibody binds human factor XII in solid-phase ELISA and immunoblotting techniques. Approximate working concentration is 1-5 ug/mL. Researcher should titer antibody in specific assay. Reconstitute with deionized water, aliquot and store at -20°C for prolonged periods. Avoid freeze-thaw cycles.

For Research Use Only.