

Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK
Tel: +44 (0)1792 686070
Email: info@enzymeresearch.co.uk
Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	Mab HPC
Lot:	SAMPLE
Description:	Monoclonal antibody to Human Protein C
Concentration:	1 mg/mL (before lyophilization)
Aliquots:	1 x 0.50 mg
Volume/aliquot:	0.500 ml
Activity:	N/A
Total Protein:	0.50 mg
Buffer:	10 mM Sodium Phosphate / 0.15 M NaCl / 0.01% Sodium Azide / pH 7.4
Storage:	-20°C
Expiration Date:	**/****

Epitope on Protein C heavy chain.

Purified from cell supernatant by Protein G affinity chromatography. Binds Protein C in solid-phase ELISA. Blots non-reduced and heavy chain of reduced human Protein C.

Reconstitution/Storage: Store lyophilized antibody at -20oC. Reconstitute with distilled water. Centrifuge to remove insoluble material. Store reconstituted antibody at 4oC or aliquot and store at -20oC. Avoid repeated freeze-thaw cycles.

For Research Use Only.