



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK
Tel: +44 (0)1792 686070
Email: info@enzymeresearch.co.uk
Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	Mab HFV
Lot:	SAMPLE
Description:	Monoclonal antibody to Human Factor V
Concentration:	1 mg/ml (before lyophilization)
Aliquots:	1 x 0.50 mg
Volume/aliquot:	0.500 mL
Activity:	N/A
Total Protein:	0.50 mg
Buffer:	10 mM Sodium Phosphate/ 0.15 M NaCl/ 0.01% Sodium Azide/ pH 7.4
Storage:	-20°C
Expiration Date:	**/****

Epitope on Factor V heavy chain.

Antibody purified by Protein G affinity chromatography from cell supernatant. Binds Human factor V in solid phase ELISA. Blots factor V and factor Va heavy chain.

We suggest hydrating the antibody to its original concentration of 1 mg/mL and aliquot in to single use size (avoid freeze/thaw cycles). Antibody expires 2 years after reconstitution.

For Research Use Only.