

Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK

Tel: +44 (0)1792 686070

Email: info@enzymeresearch.co.uk

Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	RBIIa
Lot:	SAMPLE
Description:	Rabbit Thrombin
Concentration:	0.97 mg/mL (before lyophilization)
Aliquots:	1 x 0.10 mg
Volume/aliquot:	0.103 mL
Activity:	3000.00 NIH u/mg
Total Protein:	0.10 mg
Buffer:	50 mM sodium citrate, 0.2 M NaCl, 0.1% PEG-8000, pH 6.5
Storage:	2-8 °C
Expiration Date:	**/****

The Rabbit Thrombin is purified from rabbit plasma. Purity is judged by SDS-PAGE to be >90%. We suggest hydrating the protein to the original volume with water.

For Research Use Only.