



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK

Tel: +44 (0)1792 686070

Email: info@enzymeresearch.co.uk

Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	RVV-X
Lot:	SAMPLE
Description:	Russells' Viper Venom
Concentration:	1.83 mg/mL (before lyophilization)
Aliquots:	1 x 0.10 mg
Volume/aliquot:	0.055 mL
Activity:	n/a
Total Protein:	0.10 mg (1 x 0.10 mg)
Buffer:	20 mM Tris-HCl / 0.1 M NaCl / 1 mM Benzamidine/ pH 7.4
Storage:	<-60 °C
Expiration Date:	**/****

The Russells' Viper Venom Factor X activator is >95% pure as judged by SDS-PAGE and has a high specific clotting activity.

For Research Use Only.